

## Instructions

# StarFish™

## Water & Gas/Vacuum Distribution Manifolds



Your Local Distributor

# Index

<b>Page 1</b>	<b>Index</b>
<b>Page 2</b>	<b>Introduction Warranty</b>
<b>Page 3</b>	<b>Safety Information Important WARNINGS</b>
<b>Page 4</b>	<b>Products &amp; Accessories</b>
<b>Page 5</b>	<b>Product &amp; Accessories - continued</b>
<b>Page 6</b>	<b>Set-up &amp; Operation - Water Distribution Manifolds</b>
<b>Page 7</b>	<b>Set-up &amp; Operation - Gas/Vacuum Distribution Manifold</b>
<b>Page 8</b>	<b>WARRANTY - EMAIL BACK</b>

**Thank you for purchasing your StarFish Distribution Manifold**

**Please read this Instruction Manual thoroughly before operating your unit.**

## **Introduction**

**The StarFish Distribution Manifolds form part of the StarFish family of products that enable you to perform productive heating and stirring experiments using existing glassware and stirring hotplates.**

The versatile Distribution Manifolds allow a coolant or a gas / vacuum to be evenly distributed to up to five vessels simultaneously.

### **Space saving.....**

Innovative compact design saves space.

### **Safe working.....**

Separate and easily identifiable Water and Gas / Vacuum Manifolds to prevent incorrect connection of fluid and gas supplies

### **Productivity.....**

Quick release connectors on inlet and outlets for quick and easy connection to vessels

### **Versatility.....**

All connectors feature leak proof shut off valves that automatically close when the couplings are parted, allowing manifolds to be used with any number of condensers without daisy chaining.



## **Warranty**

The Water and Gas Vacuum Manifolds include one year full parts and labour warranty from date of original purchase. Warranty will only be valid if a completed **Warranty Email Back** is returned within 1 month of date of purchase (see last page).

In the event of product failure please contact your local distributor.  
Please do not return any goods without prior agreement.

## Safety Information

The following symbols are intended to assist the user in the safe and efficient operation of the StarFish.

	<b>Warning</b> Applies when there is a possibility of personal injury.
	<b>Important Note</b> Alerts the user to important facts.

## Important WARNINGS

**Please read these instructions completely before using your Distribution Manifold.**

**Operate only in a fume cupboard with protective safety sash.**

**To avoid the build up of lime scale in the Water Distribution Manifold, please avoid the use of hard water.**

**The Manifolds are only suitable for the distribution of water (RR95500) or inert gas/vacuum (RR95510).**

**The Distribution Manifold should only be operated by trained and competent personnel.**

**As with all chemistry, care should be taken to monitor your experiment at all stages.**

**The Manifold should not be left unattended unless in a supervised area.**



### Important Note


Water Manifold RR95500 is fitted with a female inlet/outlet connector, whilst the Gas/Vacuum Manifold RR95510 is fitted with a male inlet/outlet connector. This prevents incorrect connection of fluid and gas supplies.



### Warning

The manifold can be used for both vacuum and gas purging.

The maximum operating pressure is 3psi above atmospheric pressure and a vacuum of approximately 150 to 125mBar.



### Important Note

All of the connectors feature leak-proof shut-off valves that are automatically closed when the couplings are parted. This enables the manifold to be used with any number of vessels (up to five) without daisy-chaining.



### Warning

Acetal Quick release fittings may be adversely affected by direct exposure to chemical solutions or vapours.

## Products & Accessories

Cat No.	Description	Qty
<b>StarFish Distribution Manifolds</b>		
RR95500	Water Manifold with connectors	1
RR95510	Gas/Vacuum Manifold with connectors	1
RR95520	R/Angle 6.4mm Connector + shut-off	5
RR95525	Water Manifold 9.6mm Connector + shut-off	1
RR95530	Gas Manifold 9.6mm Connector + shut-off	1
RR95535	Tubing for Manifold outlets 15m x 6.4mm	1
RR95540	Tubing for Manifold inlets 15m x 8mm	1
<b>Stirring Hotplates</b>		
RR98072	Carousel Stirring Hotplate 230v UK Plug (similar to IKA RCT)	1
RR98073	Carousel Temperature Controller (similar to IKA RCT)	1
RR95810	Cimarec 2 Digital Stir Hotplate 240v UK	1
RR95820	Heidolph MR3001K Stir Hotplate 240v UK	1
RR95825	Heidolph EKT3001 Temp Control 240v UK	1
RR95830	Stuart SB162 Stir Hotplate 240v UK	1
RR95835	Stuart SD162 Digital Stir Hotplate 240v UK	1
<b>StarFish Bases</b>		
RR95010	StarFish Base Plate 135mm (for RDT & IKA)	1
RR95020	StarFish Base Plate 145mm (for Heidolph)	1
RR95030	StarFish Base Plate 150mm (for Eyela)	1
RR95040	StarFish Base Plate 184mm x 184mm (for Cimarec)	1
RR95050	StarFish Base Plate 160mm x 160mm (for Stuart)	1
<b>Base Plate Handles</b>		
RR95100	StarFish Base Plate Handles (Pair)	1
<b>MonoBlocks</b>		
RR95130	<b>MonoBlock</b> for 5 x 250ml Flasks	1
RR95135	<b>MonoBlock</b> for 16 x 25mm Tubes	1
RR95140	<b>MonoBlock</b> for 16 x 24mm Tubes	1
RR95145	<b>MonoBlock</b> for 40 x 16mm Tubes	1
RR95150	<b>MonoBlock</b> for 40 x 12mm Tubes	1
RR95155	<b>MonoBlock</b> for 20 x 21mm Vials (4 dram)	1
RR95160	<b>MonoBlock</b> for 40 x 17mm Vials (2 dram)	1
RR95165	<b>MonoBlock</b> for 40 x 15mm Vials (1 dram)	1
RR95170	<b>MonoBlock</b> for 40 x 12mm vial (2ml)	1
<b>PolyBlocks</b>		
RR95230	<b>PolyBlock</b> for 1 x 250ml Flask	1
RR95235	<b>PolyBlock</b> for 3 x 25mm Tubes	1
RR95240	<b>PolyBlock</b> for 3 x 24mm Tubes	1
RR95245	<b>PolyBlock</b> for 9 x 16mm Tubes	1
RR95250	<b>PolyBlock</b> for 9 x 12mm Tubes	1
RR95255	<b>PolyBlock</b> for 3 x 21mm Vials (4 dram)	1
RR95260	<b>PolyBlock</b> for 7 x 17mm x Vials (2 dram)	1
RR95265	<b>PolyBlock</b> for 9 x 15mm x Vials (1 dram)	1
RR95270	<b>PolyBlock</b> for 9 x 12mm x Vials (2ml)	1
<b>StarFish Inserts</b>		
RR95330	150ml Flask Insert	1
RR95335	100ml Flask Insert	1
RR95340	50ml Flask Insert	1
RR95345	25ml Flask Insert	1
RR95350	10ml Flask Insert	1
RR95355	5ml Flask Insert	1
<b>StarFish Support Rods</b>		
RR95665	650mm Rod	1
RR95666	650mm Split Rod	1
<b>StarFish Clamp</b>		
RR95400	Universal Telescopic 5-way Clamp (with Velcro)	1
RR95430	Replacement Self Adhesive Velcro Pads	10
RR95440	Replacement Velcro Loop Strips 200mm	5
RR99613	Pivot Ring Stirring Bar 15mm	40
RR95925	Spinning Stirring Bar 31mm	1
RR95926	Spinning Stirring Bar 44mm	1
RR95927	Spinning Stirring Bar 57mm	1

## Products & Accessories

<b>Test Tubes</b>		
RR94005	24x150 T/Tube with Rim (Medium Wall)	100
RR94010	24x150 T/Tube without Rim (Medium Wall)	100
RR94015	16x100 T/Tube with Rim (Medium Wall)	100
RR94020	16x100 T/Tube without Rim (Medium Wall)	100
RR94025	12x100 T/Tube with Rim (Medium Wall)	100
RR94030	12x100 T/Tube without Rim (Medium Wall)	100
<b>Round Bottom Flasks</b>		
RR94050	250ml Round Bottom Flask 24/29 Socket	5
RR94055	250ml Round Bottom Flask 19/26 Socket	5
RR94060	150ml Round Bottom Flask 24/29 Socket	5
RR94065	150ml Round Bottom Flask 19/26 Socket	5
RR94070	100ml Round Bottom Flask 24/29 Socket	5
RR94075	100ml Round Bottom Flask 19/26 Socket	5
RR94080	50ml Round Bottom Flask 24/29 Socket	5
RR94085	50ml Round Bottom Flask 19/26 Socket	5
RR94090	25ml Round Bottom Flask 24/29 Socket	5
RR94095	25ml Round Bottom Flask 19/26 Socket	5
RR94100	10ml Round Bottom Flask 14/23 Socket	5
RR94105	10ml Round Bottom Flask 10/19 Socket	5
RR94110	5ml Round Bottom Flask 14/23 Socket	5
RR94115	5ml Round Bottom Flask 10/19 Socket	5
<b>Round Bottom Flasks with Sidearm</b>		
RR94150	250ml R/B Flask 24/29 with 19/26 Sidearm	2
RR94155	100ml R/B Flask 24/29 with 19/26 Sidearm	2
RR94160	50ml R/B Flask 14/23 with 14/23 Sidearm	2
RR94165	25ml R/B Flask 14/23 with 14/23 Sidearm	2
<b>Condensers</b>		
RR94180	Liebig Condenser 208mm 19/26	1
RR94185	Liebig Condenser 250mm 24/29	1
RR94195	Coil Condenser 165mm 19/26	1
RR94200	Coil Condenser 165mm 24/29	1
RR94210	Allihn Condenser 250mm 19/26	1
RR94215	Allihn Condenser 250mm 24/29	1
<b>Soxhlet</b>		
RR94230	20ml Soxhlet 19/26 Cone 24/29 Socket	1
RR94235	40ml Soxhlet 24/29 Cone 29/32 Socket	1
RR94240	60ml Soxhlet 24/29 Cone 34/25 Socket	1
RR94245	100ml Soxhlet 24/29 Cone 40/38 Socket	1
<b>Cat No.</b>	<b>Description</b>	<b>Qty</b>
<b>Right Angle Adapters</b>		
RR94260	Adapter B19 Cone to SQ13 Screwthread	5
RR94265	Adapter B24 Cone to SQ13 Screwthread	5
RR94270	Adapter B29 Cone to SQ13 Screwthread	5
RR94275	Adapter B34 Cone to SQ13 Screwthread	5
<b>Beakers</b>		
RR94350	250ml Beaker	10
RR94355	150ml Beaker	10
RR94360	100ml Beaker	10
RR94365	50ml Beaker	10
<b>Conical Flasks</b>		
RR94380	250ml Conical Flask Narrow Neck	10
RR94385	150ml Conical Flask Narrow Neck	10
RR94390	100ml Conical Flask Narrow Neck	10
RR94395	50ml Conical Flask Narrow Neck	10
<b>Dropping Funnels</b>		
RR94450	250ml Dropping Funnel 19/26	1
RR94455	100ml Dropping Funnel 19/26	1
RR94460	50ml Dropping Funnel 14/23	1
<b>Vials</b>		
RR94300	4 Dram Vial 21x70mm	200
RR94310	2 Dram Vial 17x60mm	200
RR94320	1 Dram Vial 15x45mm	100
RR94330	2ml Vial 12x35mm	1000

# Set-Up & Operation - Water Distribution Manifolds

## Water Distribution Manifold

Water Distribution Manifolds have been designed to allow coolant from a single source to be evenly distributed to up to five condensers and then the flow re-combined to one outlet pipe. Two manifolds are used in each system, one to distribute water to the condensers and one to collect coolant for recirculation or to drain.

1. Confirm that the manifolds are Water Manifold RR95500 rather than Gas/Vacuum Manifolds. The Water Manifold is identified by the label and fitted with a female connector on the side of the manifold body.

Important, Important,  
Important, Important,  
Important, Important,  
Important, Important,  
Important, Important,  
Important, Important,  
Important, Important,

### Important Note

Water Manifold RR95500 is fitted with a female inlet/outlet connector, whilst the Gas/Vacuum Manifold RR95510 is fitted with a male inlet/outlet connector. This prevents incorrect connection of fluid and gas supplies.

2. Place both Water Manifolds over the top of support rod and slide down to the desired position and fix with the thumbwheel. (It does not matter which way up the manifolds are fitted to the support rod). The Manifold will accept support rods up to a maximum of 12mm OD.
3. Push the the large barbed connector RR95525 (supplied with the Water Manifold) into your water inlet tubing. (For suitable tubing please order RR95540 - 8mm x 15m clear Tygon Tubing).

*NB. The RR95525 connector has a 9.6mm OD barb with an 6.4mm bore and accepts flexible tubing with an 8mm ID.*

4. Now the tubing is connected to the RR95525 connector, insert it into the female connector on the side of the lower of the two Water Manifolds. It will click in place.
5. Repeat the procedure with the water outlet/drain tubing. Connecting it to the RR95525 connector on the upper of the two Water Manifolds.
6. Connect the tubing from each of your condenser inlets to a right-angled connector RR95520 (5 are supplied with each Water Manifold). (For suitable tubing please order RR95535 - 6.4mm x 15m clear Tygon Tubing).

*NB. The RR95520 right-angled connector has a 6.4mm OD barb with an 3.2mm bore and accepts flexible tubing with an 6.4mm ID.*

7. Now the tubing is connected to the RR95520 right angled connector, insert it into any one of the five female connectors on the top of the lower Water Manifold. It will click in place.
8. Repeat the procedure with the with the upper Water Manifold, connecting each of the condenser outlets to the manifold. For best results coolant should always enter the bottom of you condenser and exit at the top.
9. A water flow rate of at least 1.5 litres per minute should be used to maintain even distribution to each condenser.

Important, Important,  
Important, Important,  
Important, Important,  
Important, Important,  
Important, Important,  
Important, Important,  
Important, Important,


### Important Note

All of the connectors feature leak-proof shut-off valves that are automatically closed when the couplings are parted. This enables the manifold to be used with any number of condensers (up to five) without daisy-chaining.

Important, Important,  
Important, Important,  
Important, Important,  
Important, Important,  
Important, Important,  
Important, Important,  
Important, Important,

### Important Note

To avoid the build-up of lime scale in the Water Distribution Manifold, please avoid the use of hard water.



# Set-Up & Operation - Gas/Vacuum Distribution Manifold

## Gas/Vacuum Distribution Manifold

Gas/Vacuum Distribution Manifolds have been designed to allow gas or a vacuum from a single source to be evenly distributed to up to five vessels.

1. Confirm that the manifold is Gas/Vacuum Manifold RR95510 rather than a Water/Coolant Manifold. The Gas/Vacuum Manifold is identified by the label and fitted with a male connector on the side of the manifold body.

Important, Important,  
Important, Important,  
Important, Important,  
Important, Important,  
Important, Important,  
Important, Important,  
Important, Important,

### Important Note

Water Manifold RR95500 is fitted with a female inlet/outlet connector, whilst the Gas/Vacuum Manifold RR95510 is fitted with a male inlet/outlet connector. This prevents incorrect connection of fluid and gas supplies.

2. Place the Gas/Vacuum Manifold over the top of support rod and slide down to the desired position and fix with the thumbwheel. (It does not matter which way up the manifold is fitted to the support rod). The Manifold will accept support rods up to a maximum of 12mm OD.
3. Push the the large barbed connector RR95530 (supplied with the Gas/Vacuum Manifold) into your gas inlet tubing. (For suitable tubing please order RR95540 - 8mm x 15m clear Tygon Tubing).

*NB. The RR95530 connector has a 9.6mm OD barb with an 6.4mm bore and accepts flexible tubing with an 8mm ID.*

4. Now the tubing is connected to the RR95525 connector, insert it into the male connector on the side of Gas/Vacuum Manifold. It will click in place.
5. Connect the tubing from each of your vessels to a right-angled connector RR95520 (5 are supplied with each Gas/Vacuum Manifold). (For suitable tubing please order RR95535 - 6.4mm x 15m clear Tygon Tubing).

*NB. The RR95520 right-angled connector has a 6.4mm OD barb with an 3.2mm bore and accepts flexible tubing with an 6.4mm ID.*

7. Now the tubing is connected to the RR95520 right angled connector, insert the connector into any one of the five female connectors on the top of the Gas/Vacuum Manifold. It will click in place.



### Warning


The manifold can be used for both vacuum and gas purging.

The maximum operating pressure is 3psi above atmospheric pressure and a vacuum of approximately 150 to 125mBar.

Important, Important,  
Important, Important,  
Important, Important,  
Important, Important,  
Important, Important,  
Important, Important,

### Important Note

All of the connectors feature leak-proof shut-off valves that are automatically closed when the couplings are parted. This enables the manifold to be used with any number of vessels (up to five) without daisy-chaining.





# Warranty – Email Back

sales@radleys.co.uk

To qualify for your warranty please complete, scan and email this form to Radleys

**Date of Purchase** .....

**Supplier's Name and Address**  
.....

**Product Batch/Serial No. (if shown)** .....

**Your Details**

**Mr Mrs Miss Ms Dr Prof**

**Name** .....

**Position** .....

**Dept** ..... **Building** .....

**Organisation** .....

**Address 1** .....

**Address 2** .....

**Town/City** ..... **County/State** .....

**Country** ..... **Post/Zip Code** .....

**Telephone** ..... **Ext** ..... **Fax** .....

**Email** ..... **Website** .....

**Type of Organisation; please tick all boxes relevant**

Academic Institution     Consumer Goods     Defence/Military/Forensic     Government     Manufacturing/Industrial     Polymers/Plastics  
 Animal Health/Zoology     Contract Lab     Environmental/Water     Hospital/Pharmacy     Nuclear/Gas/Electric     Process Engineering  
 Agrochemical     Contract Synthesis     Flavours/Fragrances     Instrum/Elect & Mech     Petrochemical/Oil     Research Institute  
 Chemical Manufacture     Cosmetics     Food/Beverages     Lab Equip Dealer/Mnf     Pharma/Biotech/API     Other.....

**Areas of Interest; please tick all boxes relevant**

Analytical Chemistry     Chromatography     Estate & Facilities     Health & Safety     Organic Chemistry     QC/QA  
 Automation/HTS     Clinical/Medical/Pathology     Food & Agriculture     Inorganic/Metallurgy     Parallel Chem/Combi-Chem     Sales & Marketing  
 Biochemistry     Construction     Formulation     Liquid Handling/MicroPlates     Polymers & Oils     Separation/SPE  
 Biological Sciences     Drug Discovery     Geology     Material Science     Process Dev/Scale-up     Support/Engineering  
 Catalysis     Environmental Health     Glassblower     Medical Devices     Process Safety/Calorimetry     Temperature Control  
 Other.....     Medicinal Chemistry     Purchasing/Stores     Veterinary

To request specific product information from Radleys please fill in below

<p><b>Benchtop and Hotplates</b></p> <p> <input type="checkbox"/> Findenser Air Condenser  <input type="checkbox"/> Heat-On Block System  <input type="checkbox"/> Cool-It Insulated Bowls  <input type="checkbox"/> StarFish Work Station  <input type="checkbox"/> Carousel Stirring Hotplates  <input type="checkbox"/> Overhead Stirrers                 </p> <p><b>Jacketed Lab Reactors</b></p> <p> <input type="checkbox"/> Reactor-Ready Lab Reactor  <input type="checkbox"/> Reactor-Ready Duo Lab Reactor  <input type="checkbox"/> Reactor-Ready Pilot Lab Reactor  <input type="checkbox"/> Custom Jacketed Reaction Systems                 </p>	<p><b>Parallel Reaction Stations</b></p> <p> <input type="checkbox"/> Carousel 12 Plus Reaction Station  <input type="checkbox"/> Cooled Carousel 12 Reaction Station  <input type="checkbox"/> Carousel 6 Plus Reaction System  <input type="checkbox"/> Cooled Carousel 6 Plus Reaction Station  <input type="checkbox"/> Carousel Work-Up Station  <input type="checkbox"/> GreenHouse Plus Parallel Synthesiser  <input type="checkbox"/> GreenHouse Work-Up Station  <input type="checkbox"/> GreenHouse Blowdown Evaporator  <input type="checkbox"/> Tornado Overhead Stirring System  <input type="checkbox"/> Breeze Heating/Cooling Work Station  <input type="checkbox"/> Storm Heating/Cooling Work Station                 </p>	<p><b>Software</b></p> <p> <input type="checkbox"/> AVA Lab Control Software  <input type="checkbox"/> Level 1/2  <input type="checkbox"/> Level 3/4  <input type="checkbox"/> Data Hub                 </p> <p><b>Automated Reaction Stations</b></p> <p> <input type="checkbox"/> Mya 4 Reaction Station                 </p> <p><b>Other</b></p> <p> <input type="checkbox"/> Huber.....  <input type="checkbox"/> Heidolph.....  <input type="checkbox"/> Other.....                 </p>
--	--	---