Biosan SIA., Ratsupites 7, build. 2, Riga, LV-1067, Latvia, Phone: +371 674 261 37, Fax: +371 674 281 01, E-mail: marketing@biosan.lv, www.biosan.lv

# HSP-30/100, platform

## DESCRIPTION

Platform with tight fit clamps for flasks.

## SPECIFICATIONS

Working area	360 × 400 mm
Flask volume	100-150 ml
Number of clamps	30



#### CAT. NUMBER

BS-010167-KK

Platform HSP-30/100







ES-20/80 Software included, without platform Orbital Shaker-Incubator

**ES-20/80** shaker-incubator for biotechnological and pharmaceutical laboratories is a professional category equipment. The typical applications include - microbial and cell culture cultivation, protein expression, solubility studies, general mixing, as well ...

#### read more



ES-20/80C Software included, without platform Orbital Shaker-Incubator with cooling

**ES-20/80C** shaker-incubator for biotechnological and pharmaceutical laboratories is a professional category equipment. The typical applications include - microbial and cell culture cultivation, protein expression, solubility studies, general mixing, as well ...

read more